


CARE GUIDELINES

- The DCI Lounger should be serviced by specialized service provider indicated by the vendor. Unplug the DCI Lounger before starting any service or maintenance work. Do not plug the DCI Lounger to a damaged electrical outlet.
- If the DCI Lounger is not going to be used during prolonged periods, coil the cable, cover and store it in a dry place. Do not put DCI Lounger near fire and protect it against prolonged exposure to the sun. Clean the DCI Lounger with a dry cloth, do not use solvents, petrol or alcohol. It is recommended to use wood cleaning and care products.
- Components have been designed and adapted specifically for the DCI Lounger. Do not move the DCI Lounger on uneven surfaces. The DCI Lounger should be moved on its wheels.
- Ambient temperature appropriate to turn the device on: -10°C to $+50^{\circ}\text{C}$ (14°F to 122°F).

COMFORT LINE



PREMIUM LINE

